

MOUNTING SELECTION SHEET

CLIENT: *The client is the person who will be receiving the equipment or services. (NOT FOR RENTAL. PURCHASE ONLY FORM)*

Client Last Name: _____ First Name: _____

Make/Model of Wheelchair: _____ Device/Access Method: _____

CONTACT PERSON: *Contact person is the person who would address any questions related to this Selection Sheet.*

Contact Last Name: _____ First Name: _____

Contact Information (phone/email): _____

All PRC-Salttillo stocked mounts are manufactured by Rehadapt. If you are requesting a mount **NOT** on this form (either from Rehadapt or another manufacturer), please attach a quote for that mount. Note: custom ordered mounts **cannot** be returned or exchanged.

Mounting Plate (E2599)



The Mounting Plate is included with the Accent® 1400 and Accent 1000 with Look®, but NOT included with other PRC-Salttillo device options. If ordering for a device without a plate and you would like a plate sent for that device, list device model and serial number (if known): _____

Step A

if only selecting wheelchair mounts, skip Step A and continue to Step B

Floorstands (E2512)



Pentalock Tele
includes counterweight



TeleLock
add quickpack option*



VarioLock
add floating option*

Table stands & clamps (E2512)



TS-XL



TS-GO



TS-UC



ClampOnMount

All Purpose Mounts (E2512) are available to accommodate variable or complex positioning of your speech-generating device in multiple locations such as wheelchair and bed/chair/floor/etc. The VarioLock with floating option and the TeleLock with quickpack option are not available for the All Purpose Mount. For pricing, configuration, and availability, contact your PRC-Salttillo consultant or the Funding Department.

**Due to the cost of this product, it may not be available for all funding sources. Contact the Funding Department with questions.*

MOUNTING SELECTION SHEET

CLIENT: *The client is the person who will be receiving the equipment or services.*

Client Last Name: _____ First Name: _____

Contact Last Name: _____ First Name: _____

Step B

continue to Step C if selecting wheelchair mounts

Wheelchair Mounts and Tubing Size (E2512)

F = frame tube, **V** = vertical tube, **H** = horizontal tube



RM-2

two tubes

Size options in inches
Vertical: 12, 16, 20, 24
Horizontal: 8, 12, 16

V	H



RM-3

three tubes

Size options in inches
Frame: 8, 12
Vertical: 12, 16, 20, 24
Horizontal: 8, 12, 16

F	V	H



RM-3-Hybrid**

thinner end tubes

Size options in inches
Frame: 8, 12
Vertical: 12, 16
Horizontal: 8, 12, 16

F	V	H

****Use RM3-Hybrid with Accent 800, NovaChat® 8, and Via® Mini. Hybrids are light weight mounts. Consider RM-2 and RM-3 for greater support and stability.**



RM-3-HD

double base tube

Size options in inches
Vertical: 12, 16, 20
Horizontal: 8, 12, 16

V	H



RM-S

curved main tube



Curved tube - choose one

Small (offset 6", height 26")

Medium (offset 9", height 30")

Large (offset 12", height 34")

Horizontal tube (H) - choose one

Ship standard tube length (16")

Replace horizontal tube (12")

Larger offset will set the device further from the client.

Wheelchair mounts are complex. We can help! If you have questions about a specific mount and/or frame clamp, please email mounting@prc-salttillo.com with front, back, left, and right side pictures of your wheelchair and they can assist you.

MOUNTING SELECTION SHEET

CLIENT: *The client is the person who will be receiving the equipment or services.*

Client Last Name: _____ First Name: _____

Contact Last Name: _____ First Name: _____


Step C

Frame clamps

Select one clamp, regardless of wheelchair mount.


Mount Location: left side • right side

for round tubes



round	size	side
round	7/8"	--
round	1"	side
round	1 1/8"	--
round	1 1/4"	--
round	1 3/8"	--

for rectangular tubes



large	
small	

for T-nut rails



Quantum/Invacare	
Quickie	
Permobil	

for frame with holes



M6/M8 Thread (~1/4" and 5/16")	
-----------------------------------	--


Additional frame clamp required? Please specify from the list above. For a quick and easy switch between two wheelchairs, we recommend the Base Coupler Option (see below).

Other frame clamps If the frame clamp you would like is not listed, please note the frame clamp name/info along with the make/model of your chair and email our mounting team at mounting@prc-salttillo.com.

Step D optional

Additional options available at an additional cost (E2512)

May not be available for all funding sources.



Base Coupler Option

Coupler with a quick release mechanism to remove the entire mount directly at the wheelchair frame



Second Coupler

Add a Second Coupler half for second wheelchair

Fold Down Joint Option for quick folding of mount

left side • right side

Left and right is from perspective of person in wheelchair. Fold Down NOT available for hybrids.



Universal Quick Shift Joint

Included if Fold Down Joint option is not selected above

Offset Spacers for additional clearance when mounting around obstacles

1" Spacer QTY: _____

1 1/2" Spacer QTY: _____

