

MOUNTING SELECTION SHEET

CLIENT: *The client is the person who will be receiving the equipment or services. (NOT FOR RENTAL. PURCHASE ONLY FORM)*

Client Last Name: _____ First Name: _____

Make/Model of Wheelchair: _____ Device/Access Method: _____

CONTACT PERSON: *Contact person is the person who would address any questions related to this Selection Sheet.*

Contact Last Name: _____ First Name: _____

Contact Information (phone/email): _____

All PRC-Salttillo stocked mounts are manufactured by Rehadapt. If you are requesting a mount **NOT** on this form (either from Rehadapt or another manufacturer), please attach a quote for that mount. Note: custom ordered mounts **cannot** be returned or exchanged.



A **mount plate** is included with the Accent 1400®, Accent 1000 with Look®, and standard mount orders for a specified device. Select the device you will be using with this mount so we can be sure to send you the correct mount plate. If your device is not listed or a mount plate is not needed, please select "None of the above" or contact your consultant for more options. *Custom mount orders do not include a mount plate, but it may be purchased separately.*

Accent 800
Via Mini

Accent 1000 (without Look)
Via Pro

NovaChat 8
None of the above

NovaChat 10

Step A

if only selecting wheelchair mounts, skip Step A and continue to Step B

Floorstands (E2512)



Pentalock Tele
includes counterweight



TeleLock
add quickpack option*



VarioLock
add floating option*

Table stands & clamps (E2512)



TS-XL



TS-GO



TS-UC



ClampOnMount

All Purpose Mounts (E2512) are available to accommodate variable or complex positioning of your speech-generating device in multiple locations such as wheelchair and bed/chair/floor/etc. The VarioLock with floating option and the TeleLock with quickpack option are not available for the All Purpose Mount. For pricing, configuration, and availability, contact your PRC-Salttillo consultant or the Funding Department.

**Due to the cost of this product, it may not be available for all funding sources. Contact the Funding Department with questions.*

MOUNTING SELECTION SHEET

CLIENT: *The client is the person who will be receiving the equipment or services.*

Client Last Name: _____ First Name: _____

Contact Last Name: _____ First Name: _____

Step B

continue to Step C if selecting wheelchair mounts

Wheelchair Mounts and Tubing Size (E2512)

F = frame tube, V = vertical tube, H = horizontal tube



RM-2
two tubes
Size options in inches
Vertical: 12, 16, 20, 24
Horizontal: 8, 12, 16

V	H



RM-3
three tubes
Size options in inches
Frame: 8, 12
Vertical: 12, 16, 20, 24
Horizontal: 8, 12, 16

F	V	H



RM-3-Hybrid**
thinner end tubes
Size options in inches
Frame: 8, 12
Vertical: 12, 16
Horizontal: 8, 12, 16

F	V	H

***Use RM3-Hybrid with Accent 800, NovaChat® 8, and Via® Mini. Hybrids are light weight mounts. Consider RM-2 and RM-3 for greater support and stability.*



RM-3-HD
double base tube
Size options in inches
Vertical: 12, 16, 20, 24
Horizontal: 8, 12, 16

V	H



RM-S
curved main tube



Curved tube - choose one
Small (offset 6", height 26")
Medium (offset 9", height 30")
Large (offset 12", height 34")
Horizontal tube (H) - choose one
Ship standard tube length (16")
Replace horizontal tube (12")

Larger offset will set the device further from the client.

Wheelchair mounts are complex. We can help! If you have questions about a specific mount and/or frame clamp, please email mounting@prc-salttillo.com with front, back, left, and right side pictures of your wheelchair and they can assist you.

MOUNTING SELECTION SHEET

CLIENT: *The client is the person who will be receiving the equipment or services.*

Client Last Name: _____ First Name: _____


Contact Last Name: _____ First Name: _____

Step C

Frame clamps

Select one clamp, regardless of wheelchair mount.

for round tubes



round	size	side
round	7/8"	--
round	1"	side
round	1 1/8"	--
round	1 1/4"	--
round	1 3/8"	--

for rectangular tubes



large	
small	

for T-nut rails



Quantum/Invacare	
Quickie	
Permobil	

for frame with holes




M6/M8 Thread (~1/4" and 5/16")	
--------------------------------	--

Additional frame clamp required? Please specify from the list above. For a quick and easy switch between two wheelchairs, we recommend the Base Coupler Option (see below).


Other frame clamps If the frame clamp you would like is not listed, please note the frame clamp name/info along with the make/model of your chair and email our mounting team at mounting@prc-salttillo.com.

Step D optional


Additional options available at an additional cost (E2512) May not be available for all funding sources.




Base Coupler Option
 Coupler with a quick release mechanism to remove the entire mount directly at the wheelchair frame




Second Coupler
 Add a Second Coupler half for second wheelchair



Fold Down Joint Option for quick folding of mount
 left side • right side
 Left and right is from perspective of person in wheelchair.
 Fold Down NOT available for hybrids.



Universal Quick Shift Joint
Included if Fold Down Joint option is not selected above



Offset Spacers for additional clearance when mounting around obstacles
 1" Spacer QTY: _____
 1 1/2" Spacer QTY: _____